KAPA3G Hotstart DNA Polymerase

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

45 cycles: denaturation: 1 sec at 95 °C, annealing/extension: 10 sec at 60 °C

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RT-PCR 条件: RT hold: 5 min at 55 °C, initial denaturation hold: 1 min at 95 °C,

45 cycles: denaturation: 10 s at 95 °C, annealing/extension: 30 s at 60 °C

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Ready for faster turnaround times?

Scan the QR code to see how KAPA3G HotStart DNA Polymerase and other Roche CustomBiotech enzymes can accelerate your assay. Detailed data on KAPA3G performance (speed, robustness, inhibitor tolerance) is available upon request.

Regulatory disclaimers, warranty limitations and legal notes

are listed on the respective product page on custombiotech.roche.com

Trademarks

KAPA and NXTSCRIPT are trademarks of Roche.

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Tel:0120-600-152

References

- ¹ Sexually transmitted infections (STIs). [https://www.who.int/news-room/fact-sheets/detail/sexually-transmitted-infections-\(stis\)](https://www.who.int/news-room/fact-sheets/detail/sexually-transmitted-infections-(stis)). Accessed: 10-Oct-2023
- ² Institute of Medicine (US) Committee on Prevention and Control of Sexually Transmitted Diseases. The Hidden Epidemic: Confronting Sexually Transmitted Diseases: Summary. Eng TR, Butler WT, editors. Washington (DC): National Academies Press (US); 1997. PMID: 25165803.
- ³ Sexually Transmitted Infections Prevalence, Incidence, and Cost Estimates in the United States. <https://www.cdc.gov/std/statistics/prevalence-incidence-cost-2020.htm>. Accessed: 10-Oct-2023.
- ⁴ The role of diagnostic testing in an era of antimicrobial stewardship. <https://diagnostics.roche.com/global/en/article-listing/health-topics/infectious-diseases/amr-ams.html>. Accessed 10-Oct-2023
- ⁵ PCR inhibitors – occurrence, properties and removal. Schrader C, et al. <https://ami-journals.onlinelibrary.wiley.com/doi/10.1111/j.1365-2672.2012.05384.x>. Accessed 26-Oct-2023
- ⁶ Inhibitory effects of urine on the polymerase chain reaction for cytomegalovirus DNA. <https://pubmed.ncbi.nlm.nih.gov/1646235/>. Accessed 26-Oct-2023

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All rights reserved.