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Galactosylation of IgG

1. Preparation of the reaction buffer, containing 10 mM UDP-Gal, 5 mM Manganese(II) Chloride (MnCl₂), and 100 mM MES, pH 6.5.
 - Dissolve 61 mg UDP-Gal, 6.3 mg MnCl₂, and 195 mg MES in 6 ml water.
 - Adjust pH to 6.5, and adjust to a final volume of 10 ml with water.
2. Setting up the galactosylation reaction, containing 976 µg UDP-Galactose, 10 µg GalT1, and 1 mg IgG.
 - 160 µl reaction buffer
 - 33.3 µl of IgG (30 mg/ml, pH 6)
 - 1.8 µl Galactosyltransferase
 - 4.9 µl water
3. Incubate the reaction at 37°C, then take aliquots and stop the reaction by freezing the tube at -15°C to -25°C for mass spectrometry analysis.

Sialylation of IgG using either ST6 or ST3

1. Preparation of the reaction buffer, containing 6.1 mM CMP-NANA, 200 mM MES, pH6.5.
 - Dissolve 40 mg CMP-NANA and 390.4 mg MES in 6 ml of water.
 - Adjust pH to 6.5, and adjust to a final volume of 10 ml with water.
2. Setting up the sialylation reaction, containing 500 µg CMP-NANA, 100 µg Sialyltransferase, and 1 mg IgG.
 - 125 µl reaction buffer
 - 100 µl IgG (10 mg/ml, pH 6)
 - 18 µl Sialyltransferase (ST3 or ST6)
 - 2.8 µl Alkaline Phosphatase (1 mg/ml)
 - 2.53 µl ZnCl₂ solution (10 mM)
 - 1.67 µl water
3. Incubate the reaction at 37°C, then take aliquots and stop the reaction by freezing the tube at -15°C to -25°C for mass spectrometry analysis.

Simultaneous galactosylation and sialylation (“1-pot-reaction”)

1. Setting up the galactosylation reaction, containing 976 µg UDP-Galactose, 10 µg GalT1, and 1 mg IgG.
 - 160 µl reaction buffer
 - 33.3 µl of IgG (30 mg/ml, pH 6)
 - 1.8 µl Galactosyltransferase
 - 4.9 µl water
2. Incubate the reaction for 4 hours at 37°C, then take 50 µl of this reaction for mass spectrometry analysis.
3. Add reagents for subsequent sialylation to the remaining 150 µl galactosylated IgG, to set up the sialylation reaction, containing 375 µg CMP-NANA, 75 µg Sialyltransferase, and 0.75 mg IgG.
 - 94 µl reaction buffer
 - 13.6 µl Sialyltransferase ST6
 - 2.9 µl Alkaline Phosphatase (1 mg/ml)
 - 2.6 µl ZnCl₂ solution (10 mM)
4. Incubate the reaction at 37°C, then take aliquots and stop the reaction by freezing the tube at -15°C to -25°C for analysis.

Note: In case of varying IgG concentrations, the protocol must be adapted. The most important parameter to be considered is the mass ratio of activated sugar, IgG and glycosyltransferase.

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Regulatory disclaimers and license limitations

are listed on the respective product pages on custombiotech.roche.com

References

¹ Post-translational modifications in the context of therapeutic proteins, Gary Walsh & Roy Jefferis, Nat Biotechnol. 2006, Oct;24(10):1241-52

² In vitro glycoengineering of IgG1 and its effect on Fc receptor binding and ADCC activity, M. Thomann et al., PLoS One. 2015 Aug 12;10(8):e0134949

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© 2026 All rights reserved.